


## Anti-Decubitus Air Mattresses with Pump

Good for sickbed: incapable of movement, but mattresses can move

Using the method of anti bed sore cushion on a regular basis two inflatable balloon rotation, deflation, which led to parts of the body in patients with continuous change, can play a role in the artificial massage, can promote blood circulation, prevent muscle atrophy, patients do not need regular turn nor long bed sore! And the air cushion works continuously, without manual intervention!

### Breathable Althernating Pressure



**Model: ADM-001**

Long time lie bed patient: movement by mattresses

Intended for:

- 1, elderly bedridden elderly;
- 2, a patient who is bedridden because of an accident or illness
- 3, for other reasons, resulting in less activity, bedridden people longer

There is a lot of micro air jet on the mattress. The advantage is mattress surface fluctuation, with ventilation, shifting the body, it can play a role in the prevention of bed sore by long time use



**Model: ADM-001**



**Dry and comfortable**



**Exclude moisture**



**Keep dry**



**Air circulation**

The surface of the air cushion is filled with many tiny jets, which automatically control the jet stream and keep the pressure on the part of the patient in good condition